



Recombinant Human Reticulocalbin-3 (RCN3)

Product Code	CSB-YP856914HU
Abbreviation	RCN3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q96D15
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	KPSPDAGPHG QGRVHQAAPL SDAPHDDAHG NFQYDHEAFL GREVAKEFDQ LTPEESQARL GRIVDRMDRA GDGDGWVSLA ELRAWIAHTQ QRHIRDSVSA AWDTYDTRD GRVGWEELRN ATYGHYAPGE EFHDVEDAET YKKMLARDER RFRVADQDGD SMATREELTA FLHPEEFPHM RDIVIAETLE DLDRNKDGYV QVEEYIADLY SAEPEGEEPA WVQTERQQFR DFRDLNKDGH LDGSEVGHVW LPPAQDQPLV EANHLLHESD TDKDGRLSKA EILGNWNMFV GSQATNYGED LTRHHDEL
Source	Yeast
Target Names	RCN3
Protein Names	Recommended name: Reticulocalbin-3 Alternative name(s): EF-hand calcium-binding protein RLP49
Expression Region	21-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>