



Recombinant Human RaiBP1-associated Eps domain-containing protein 1 (REPS1)

Product Code	CSB-BP856915HU
Abbreviation	REPS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96D71
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGLTLSDAE QKYYSDFSY CDIESTKKVV VNGRVLELFR AAQLPNDVVL QIMELCGATR LGYFGRSQFY IALKLVAVAQ SGFPLRVESI NTVKDLPLPR FVASKNEQES RHAASYSSDS ENQGSYSGVI PPPPGRGQVK KGSVSHDTVQ PRTSADAQEP ASPVVSPQQS PPTSPHTWRK HSRHPSSGNS ERPLAGPGPF WSPFGEAQSG SSAGDAVWSG HSPPPPQENW VSFADTPPTS TLLTMHPASV QDQTTVRTVA SATTAEIRR QSSSYDDPWK ITDEQRQYYV NQFKTIQPD L NGFIPGSAK EFFTKSKLPI LELSHIWELS DFDKDGALT DEFCAAFHLV VARKNGYDLP EKLPELMPK LIDLESDADV GDQPGEVGYS GSPAEAPPSK SPSMPSLNQT WPELNQSSEQ WETFSESSS SQTLTQFDSN IAPADPDTAI VHPVPIRMTP SKIHMQEMEL KRTGSDHTNP TSPLLVKPSD LLEENKINSS VKFASGNTVA DGYSSSDSFT SDPEQIGSNV TRQRSHSGTS PDNTAPPPPP PRPQPSHSRS SSLDMNRTFT VTTGQQQAGV VAHPPAVPPR PQPSQAPGPA VHRPVDADGL ITHTSTSPQQ IPEQPNFADF SQFEVFAASN VNDEQDDEAE KHPEVLPAEK ASDPASSLRV AKTDSKTEEK TAASAPANVS KGTTPLAPPP KPVRRRLKSE DELRPEVDEH TQKTGVLA AV LASQPSIPRS VGKDKKAIQA SIRRNKETNT VLARLNSELQ QQLKDVLEER ISLEVQLEQL RPFSHL
Source	Baculovirus
Target Names	REPS1
Protein Names	Recommended name: RaiBP1-associated Eps domain-containing protein 1 Alternative name(s): RaiBP1-interacting protein 1
Expression Region	1-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.