



Recombinant Human E3 ubiquitin-protein ligase TRIM11 (TRIM11)

Product Code	CSB-MP856931HU
Abbreviation	TRIM11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96F44
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPDLSTNL QEEATCAICL DYFTDPVMTD CGHNFCRECI RRCWGQPEGP YACPECRELS PQRNLRPNRP LAKMAEMARR LHPPSPVPQG VCPAHREPLA AFCGDELRLI CAACERSGEH WAHRVRPLQD AAEDLKAKLE KSLEHLRKQM QDALLFQAQA DETCVLWQKM VESQRQNVLG EFERLRRLLA EEEQQLLQRL EEEELEVLPR LREGAAHLGQ QSAHLAELIA ELEGRCQLPA LGLLQDIKDA LRRVQDVKLQ PPEVVPME LR TVCRVPLVE TLRRFRGDVT LDPDTANPEL ILSEDRRSVQ RGD LRQALPD SPERFDPGPC VLGQERFTSG RHYWEVEVGD RTSWALGVCR ENVNRKEKGE LSAGNGFWIL VFLGSYNNSS ERALAPLRDP PRRVGIFLDY EAGHLSFYSA TDGSLIFIP EIPFSGTLRP LFSPLSSSPT PMTICRPKGG SGGDTLAPQ
Source	Mammalian cell
Target Names	TRIM11
Protein Names	Recommended name: E3 ubiquitin-protein ligase TRIM11 EC= 6.3.2.- Alternative name(s): Protein BIA1 RING finger protein 92 Tripartite motif-containing protein 11
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the nucleus and the cytoplasm. Its function has not been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.