



Recombinant Human E3 ubiquitin-protein ligase pellino homolog 1 (PELI1)

Product Code	CSB-MP856932HU
Abbreviation	PELI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FA3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFSPDQENHP SKAPVKYGEL IVLGYNGSLP NGDRGRRKSR FALFKRPKAN GVKPVSTVHIA CTPQAAKAIS NKDQHSISYT LSRAQTVVVE YTHDSNTDMF QIGRSTESPI DFVVTDTVPG SQSNSDTQSV QSTISRFACR IICERNPPFT ARIYAAGFDS SKNIFLGEKA AKWKTS DGQM DGLTTNGVLV MHPRNGFTED SKPGIWREIS VCGNVFSLRE TRSAQQRGKM VEIETNQLQD GSLIDL CGAT LLWRTAEGLS HTPTVKHLEA LRQEINAARP QCPVGFNTLA FPSMKRKDVV DEKQPWVYLN CGHVHGYHNW GNKEERDGKD RECPMCRSVG PYVPLWLGCE AGFYVDAGPP THAFSPCGHV CSEKTTAYWS QIPLPHGTHT FHAACPFAH QLAGEQGYIR LIFQGPLD
Source	Mammalian cell
Target Names	PELI1
Protein Names	Recommended name: E3 ubiquitin-protein ligase pellino homolog 1 Short name= Pellino-1 EC= 6.3.2.- Alternative name(s): Pellino-related intracellular-signaling molecule
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.