



Recombinant Human Uncharacterized protein LOC113230

Product Code	CSB-EP856933HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FF7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPQRAGPGL GKPSSKKSRA LRGGPGVPRA RQRSPRAAAE LNTEVPKPSA KDVEILEQES EVLEPSVEQD AESPGLSTAA GFESPEPGRD DPGPSEAQRA EVQGQHPPPG PKTSSSRSGS EAPHLRFPPC PEAPQRPPEE RFMETPIERE IRRSCEREES LRRSRGLSPG RAGRELVLR VRPVLNLPGP GPALPRALER ARAGAQMQRD IEREAHRQAA LARPAVPEPR ARSPPQPLGE LKRFFEEAAG SGSSAGAGDG AGPQRLPEPG GRPRSAVQGG CRVLGSAPPP FTPSLLEQEV RAVREREQEL QRQRRSVYGT AEFKEPTPSL TASRGDGKLV VIWPPRRKVS ENGLEQEERK P
Source	E.coli
Target Names	MISP3
Protein Names	Recommended name: Uncharacterized protein LOC113230
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.