



# Recombinant Human Uncharacterized protein LOC113230

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP856933HU                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                |
| <b>Uniprot No.</b>       | Q96FF7                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | MAPQRAGPGL GKPSSKKSRA LRGGPGVPRA RQRSPRAAAE LNTEVPKPSA<br>KDVEILEQES EVLEPSVEQD AESPGLSTAA GFESPEPGRD DPGPSEAQRA<br>EVQGQHPPPG PKTSSSRSGS EAPHLRFPPC PEAPQRPPEE RFMETPIERE<br>IRRSCEREES LRRSRGLSPG RAGRELVLR VRPVLNLPGP GPALPRALER<br>ARAGAQMQRD IEREAHRQAA LARPAVPEPR ARSPPQPLGE LKRFFEEAAG<br>SGSSAGAGDG AGPQRLPEPG GRPRSAVQGG CRVLGSAPPP<br>FTPSLLEQEV RAVREREQEL QRQRRSVYGT AEFKEPTPSL TASRGDGKLV<br>VIWPPRRKVS ENGLEQEERK P |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Target Names</b>      | MISP3                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Uncharacterized protein LOC113230                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-361                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                               |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                       |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                |