



Recombinant Human Secernin-2 (SCRN2)

Product Code	CSB-EP856936HU
Abbreviation	SCRN2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q96FV2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASSSPDSPC SCDCFVSVPP ASAIPAVIFA KNSDRPRDEV QEVVFPAGT HTPGSRLQCT YIEVEQVSKT HAVILSRPSW LWGAEMGANE HGVCIGNEAV WTKEPVGEGE ALLGMDLLRL ALERSSSAQE ALHVITGLLE HYGQGGNCLE DAAPFSYHST FLLADRTEAW VLETAGRLWA AQRIQEGARN ISNQLSIGTD ISAQHPELRT HAQAKGWWDG QGAFDFAQIF SLTQQPVRME AAKARFQAGR ELLRQRQGGI TAEVMMGILR DKESGICMDS GGFRTTASMV SVLPQDPTQP CVHFLTATPD PSRSVFKPFI FGMGVAQAPQ VLSPTFGAQD PVRTLPRFQT QVDRRHLYR GHQAALGLME RDQDRGQQLQ QKQQDLEQEG LEATQGLLAG EWAPPLWELG SLFQAFVKRE SQAYA
Source	E.coli
Target Names	SCRN2
Protein Names	Recommended name: Secernin-2
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>