



Recombinant Human Protein RMD5 homolog B (RMND5B)

Product Code	CSB-MP856939HU
Abbreviation	RMND5B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96G75
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEQCACVERE LDKVLQKFLT YGQHCERSLE ELLHYVGQLR AELASAALQG TPLSATLSLV MSQCCRKIKD TVQKCLASDHK DIHSSVSRVG KAIDRNFDSE ICGVVSDAVW DAREQQQQIL QMAIVEHLYQ QGMLSVAEEL CQUESTLNVDL DFKQPFLELN RILEALHEQD LGPALEWAVS HRQRLELNS SLEFKLHRLH FIRLLAGGPA KQLEALSAR HFQPFARLHQ REIQVMMGSL VYLRLGLEKS PYCHLLDSSH WAEICETFTR DACSLLGLSV ESPLSVSFAS GCVALPVLMM IKAVIEQRQC TGVWNHKDEL PIEIELGMKC WYHSVFACPI LRQQTSDSNP PIKLICGHVI SRDALNKLIN GGKLCPCYCP MEQNPADGKR IIF
Source	Mammalian cell
Target Names	RMND5B
Protein Names	Recommended name: Protein RMD5 homolog B
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.