



Recombinant Human Vacuolar-sorting protein SNF8 (SNF8)

Product Code	CSB-MP856945HU
Abbreviation	SNF8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96H20
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MHRRGVGAGA IAKKKLAEAK YKERGTVLAE DQLAQMSKQL DMFKTNLEEF ASKHKQEIRK NPEFRVQFQD MCATIGVDPL ASGKGFWSEM LGVGDFYYEL GVQIIEVCLA LKHRNGGLIT LEELHQQVLK GRGKFAQDVS QDDLIRAIAKK LKALGTGFGI IPVGGTYLIQ SVPAELNMDH TVVLQLAEKN GYVTVSEIKA SLKWETERAR QVLEHLLKEG LAWLDLQAPG EAHYWLPALF TDLYSQEITA EEAREALP
Source	Mammalian cell
Target Names	SNF8
Protein Names	Recommended name: Vacuolar-sorting protein SNF8 Alternative name(s): ELL-associated protein of 30 kDa ESCRT-II complex subunit VPS22 Short name=hVps22
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.