



# Recombinant Human Probable asparagine--tRNA ligase, mitochondrial (NARS2)

<b>Product Code</b>	CSB-MP856950HU
<b>Abbreviation</b>	NARS2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96I59
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSAPFP KHKPSAKLSV RDALGAQNAS GERIKIQGWI RSVRSQKEVL FLHVNDGSSL ESLQVVADSG LDSRELNFGS SVEVQGQLIK SPSKRQNVEL KAEEKIVIGN CDAKDFPIKY KERHPLEYLR QYPHFRCRTN VLGSILRIRS EATAAIHSFF KDSGFVHIHT PIITSNDSEG AGELFQLEPS GKLKVPEENF FNVPAFLTVS GQLHLEVMSG AFTQVFTFGP TFRAENSQSR RHLAEFYMIE AEISFVDSLQ DLMQVIEELF KATTMMVLSK CPEDVELCHK FIAPGQKDRL EHMLKNNFLI ISYTEAVEIL KQASQNFTFT PEWGADLRTE HEKYLVKHCG NIPVFINYP LTLKPFYMRD NEDGPQHTVA AVDLLVPGVG ELFGGGLREE RYHFLEERLA RSGLTEVYQW YLDLRRFGSV PHGGFGMGFE RYLQCILGVD NIKDVIPFPR FPHSCLL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NARS2
<b>Protein Names</b>	Recommended name: Probable asparagine--tRNA ligase, mitochondrial EC=6.1.1.22 Alternative name(s): Asparaginyl-tRNA synthetase Short name= AsnRS
<b>Expression Region</b>	15-477
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.