



Recombinant Human 4-hydroxyphenylpyruvate dioxygenase-like protein (HPDL)

Product Code	CSB-BP856953HU
Abbreviation	HPDL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96IR7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPALRLCH IAFHVPAGQP LARNLQRLFG FQPLASREVD GWRQLALRSG DAVFLVNEGA GSGEPLYGLD PRHAVPSATN LCFDVADAGA ATRELAALGC SVPVPPVRVR DAQGAATYAV VSSPAGILSL TLLERAGYRG PFLPGFRPVS SAPGPGWVSR VDHLTLACTP GSSPTLLRWF HDCLGFCHLP LSPGEDPELG LEMTAGFGLG GLRLTALQAQ PGSIVPTLVL AESLPGATTR QDQVEQFLAR HKGPGLQHVG LYTPNIVEAT EGVATAGGQF LAPPGAYYQQ PGKERQIRAA GHEPHLLARQ GILLDGDKGK FLLQVFTKSL FTEDTFFLEL IQRQGATGFG QGNIRALWQS VQEQSARSQE A
Source	Baculovirus
Target Names	HPDL
Protein Names	Recommended name: 4-hydroxyphenylpyruvate dioxygenase-like protein EC=1.13.-.- Alternative name(s): Glyoxalase domain-containing protein 1
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.