



Recombinant Human Pleckstrin homology domain-containing family A member 8 (PLEKHA8)

Product Code	CSB-BP856957HU
Abbreviation	PLEKHA8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96JA3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGVLYKWTN YLSGWQPRWF LLCGGILSY Y DSPEDAWKGC KGSIQMAVCE IQVHSVDNTR MDLIIPGEQY FYLKARVAE RQRWLVALGS AKACLTDSRT QKEKEFAENT ENLTKMSEL RLYCDLLVQQ VDKTKEVTTT GVSNSEEGID VGTLLKSTCN TFLKTLEECM QIANAAFTSE LLYRTPPGSP QLAMLKSSKM KHPIPIHNS LERQMELSTC ENGLNMEIN GEEEILMKNK NSLYLKSAEI DCSISSEENT DDNITVQGEI RKEDGMENLK NHDNNLTQSG SDSSCSPECL WEEGKEVIPT FFSTMNTSFS DIELEDSEGI PTEAFLASCY AVVPVLDKLG PTVFAPVKMD LVGNIKKVNQ KYITNKEEFT TLQKIVLHEV EADVAQVRNS ATEALLWLKR GLKFLKGFLT EVKNGEKDIQ TALNNAYGKT LRQHHGWVVR GVFALALRAA PSYEDFVAAL TVKEGDHQKE AFSIGMQRD SLYLPAMEKQ LAILDTLYEV HGLESEDEVV
Source	Baculovirus
Target Names	PLEKHA8
Protein Names	Recommended name: Pleckstrin homology domain-containing family A member 8 Short name= PH domain-containing family A member 8 Alternative name(s): Phosphatidylinositol-four-phosphate adapter protein 2 Short name= FAPP-2 Short name
Expression Region	1-519
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.