



# Recombinant Human Pleckstrin homology domain-containing family A member 8 (PLEKHA8)

<b>Product Code</b>	CSB-EP856957HU
<b>Abbreviation</b>	PLEKHA8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96JA3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEGVLYKWTN YLSGWQPRWF LLCGGILSY Y DSPEDAWKGC KGSIQMAVCE IQVHSVDNTR MDLIIPGEQY FYLKARVAE RQRWLVALGS AKACLTDSRT QKEKEFAENT ENLKTMSSEL RLYCDLLVQQ VDKTKEVTTT GVSNSEEGID VGTLLKSTCN TFLKTLEECM QIANAAFTSE LLYRTPPGSP QLAMLKSSKM KHPIPIHNS LERQMELSTC ENGLNMEIN GEEEILMKNK NSLYLKSAEI DCSISSEENT DDNITVQGEI RKEDGMENLK NHDNNLTQSG SDSSCSPECL WEEGKEVIPT FFSTMNTSFS DIELLEDSCI PTEAFLASCY AVVPVLDKLG PTVFAPVKMD LVGNIKKVNQ KYITNKEEFT TLQKIVLHEV EADVAQVRNS ATEALLWLKR GLKFLKGFLT EVKNGEKDIQ TALNNAYGKT LRQHHGWVVR GVFALALRAA PSYEDFVAAL TVKEGDHQKE AFSIGMQRD SLYLPAMEKQ LAILDTLYEV HGLESEDEVV
<b>Source</b>	E.coli
<b>Target Names</b>	PLEKHA8
<b>Protein Names</b>	Recommended name: Pleckstrin homology domain-containing family A member 8 Short name= PH domain-containing family A member 8 Alternative name(s): Phosphatidylinositol-four-phosphate adapter protein 2 Short name= FAPP-2 Short name
<b>Expression Region</b>	1-519
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.