



Recombinant Human Ligand-dependent corepressor (LCOR)

Product Code	CSB-MP856960HU
Abbreviation	LCOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96JN0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQRMIQQFAA EYTSKNSSTQ DPSQPNSTKN QSLPKASPVT TSPTAATTQN PVLSKLLMAD QDSPLDLTVR KSQSEPSEQD GVLDLSTKKS PCAGSTSLSH SPGCSSTQGN GRPGRPSQYR PDGLRSGDGV PPRSLQDGTR EGFGHSTSLK VPLARSLQIS EELLSRNQLS TAASLGPSGL QNHGQHLILS REASWAKPHY EFNLSRMKFR GNGALSNI SD LPFLAENSAF PKMALQAKQD GKKDVSHSSP VDLKIPQVRG MDLSWESRTG DQYSYSSLVM GSQTESALSK KLRAILPKQS RKSMLDAGPD SWGSDAEQST SGQPYPTSDQ EGDPGSKQPR KKRGRYRQYN SEILEEAISV VMSGKMSVSK AQSIYGIPHS TLEYKVKERL GTLKNPPKKK MKLMRSEGPV VSVKIELDPQ GEAAQSANES KNE
Source	Mammalian cell
Target Names	LCOR
Protein Names	Recommended name: Ligand-dependent corepressor Short name= LCoR Alternative name(s): Mblk1-related protein 2
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.