



Recombinant Human PDZ and LIM domain protein 2 (PDLIM2)

Product Code	CSB-MP856962HU
Abbreviation	PDLIM2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96JY6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALTVDVAGP APWGFRTGG RDFHTPIMVT KVAERGKAKD ADLRPGDIIV AINGESAEGM LHAEAQSKIR QSPSPLRLQL DRSQATSPGQ TNGDSSLEVL ATRFQGSVRT YTESQSSLRS SYSSPTSLSP RAGSPFSPPP SSSSLTGEAA ISRSFQSLAC SPGLPAADRL SYSGRPGSRQ AGLGRAGDSA VLVLPSPSPG RSSRPSMDSE GGSLLLEDS EVFKMLQENR EGRAAPRQSS SFRLQEALE AEERGGTPAF LPSSLSPQSS LPASRALATP PKLHTCEKCS TSIANQAVRI QEGRYRHPGC YTCADCGLNL KMRGHFWVD ELYCEKHARQ RYSAPATLSS RA
Source	Mammalian cell
Target Names	PDLIM2
Protein Names	Recommended name: PDZ and LIM domain protein 2 Alternative name(s): PDZ-LIM protein mystique
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.