



Recombinant Human Calcineurin subunit B type 2 (PPP3R2)

Product Code	CSB-BP856978HU
Abbreviation	PPP3R2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LZ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GNEASYPAE MCSHFDNDEI KRLGRRFKKL DLDKSGSLSV EEFMSLPELR HNPLVRRVID VFDTDGDGEV DFKEFILGTS QFSVKGDEEQ KLRFAFSIYD MDKDGYSNG ELFQVLKMMV GNNLTDWQLQ QLVDKTIIL DKDGDGKISF EEFSAVVRDL EIHKKLVIV
Source	Baculovirus
Target Names	PPP3R2
Protein Names	Recommended name: Calcineurin subunit B type 2 Alternative name(s): Calcineurin B-like protein Short name= CBLP Calcineurin BII Short name= CNBII PPP3R1-like Protein phosphatase 2B regulatory subunit 2 Protein phospho
Expression Region	2-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.