



Recombinant Human Coiled-coil domain-containing protein 114 (CCDC114)

Product Code	CSB-BP856980HU
Abbreviation	CCDC114
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M63
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGERRAYSK EVHQRINKQL EEIRRLEEV R GDLQVQISAA QNQVKRLRDS QRLENMDRLL KGRAQVQAEI EELQEQTRAL DKQIQEWETR IFTHSKNVRS PGFILDQKVK IRRRIRILEN QLDRVTCHFD NQLVRNAALR EELDLLRIDR NRYLNVDRKL KKEIHHLHHL VSTLILSSTS AYAVREEAKA KMGLLRERAE KEEAQSEMEA QVLQRQILHL EQLHHFLKLNDRQPDVDV LEKREKQAGE VAEGVWKTSQ ERLVLCYEDA LNKLSQLMGE SDPDLLVQKY LEIEERNFAE FNFINEQNLE LEHVQEEIKE MQEALVSARA SKDDQHLLQE QQQKVLQQRM DKVHSEAERL EARFQDVRGQ LEKLNKADIQL LFTKAHCDSS MIDDLLGVKT SMGDRDMGLF LSLIEKRLVE LLTVQAFLHA QSFTSLADAA LLVLGQSLED LPKKMAPLQP PDTLEDPPGF EASDDYPMSR EELLSQVEKL VELQEQAFAQ RQKDLAAAAA KLDGTLSDVDL ASTQRAGSST VLVPTRHPPH IPGSILSHKT SRDRGSLGHV TFGGLSSSTG HLP SHITHGD PNTGHVTFGS TSASSGGHVT FRPVSASSYL GSTGYVGSSR GGENTEGGVE SGGTASDSSG GLGSSRDHVS STGPASSTGP GSSTSKDSRG
Source	Baculovirus
Target Names	CCDC114
Protein Names	Recommended name: Coiled-coil domain-containing protein 114
Expression Region	1-670
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.