



Recombinant Human Coiled-coil domain-containing protein 11 (CCDC11)

Product Code	CSB-BP856981HU
Abbreviation	CCDC11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M91
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MYSQRFGTVQ REVKGPTPKV VIVRSKPPKG QGAEHHLERI RRSHQKHNAI LASIKSSERD RLKAEWDQHN DCKILDSLVR ARIKDAVQGF IINIEERRNK LRELLALEEN EYFTEMQLKK ETIEEKKDRM REKTKLLKEK NEKERQDFVA EKLDQQFRER CEELRVELLS IHQKKVCEER KAQIAFNEEL SRQKLVEEQM FSKLWEEDRL AKEKREAQEA RRQKELMENT RLGLNAQITS IKAQRQATQL LKEEEARLVE SNNAQIKHEN EQDMLKKQKA KQETRILQK ALQERIEHIQ QEYRDEQDLN MKLVQRALQD LQEEADKKKQ KREDMIREQK IYHKYLAQRR EEKAQEKEF DRILEEDKAK KLAEKDKELR LEKEARRQLV DEVMCTRKLQ VQEKLQREAK EQEERAMEQK HINESLKELN CEEKENFARR QRLAQEYRKQ LQMQUIAYQQQ SQEAEKEEKR REFEAGVAAN KMCLDKVQEV LSTHQVLPQN IHPMRKACPS KLPP
Source	Baculovirus
Target Names	CFAP53
Protein Names	Recommended name: Coiled-coil domain-containing protein 11
Expression Region	1-514
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.