



Recombinant Human Zinc finger protein 42 homolog (ZFP42)

Product Code	CSB-MP856985HU
Abbreviation	ZFP42
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MM3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSQQLKKRAK TRHQKGLGGR APSGAKPRQG KSSQDLQAEI EPVSAVWALC DGYVCYEPGP QALGGDDFSD CYIECVIRGE FSQPILEEDS LFESLEYLKK GSEQQLSQKV FEASSLECSL EYMKKGVKKE LPQKIVGENS LEYSEYMTGK KLPPGGIPGI DLSDPKQLAE FARKKPPINK EYDSLAIAC PQSGCTRKLR NRAALRKHLL IHGPRDHVCA ECGKAFVSS KLKRHFLVHT GEKPFRCFFE GCGKRFSLDF NLRTHVRIHT GEKRFVCPFQ GCNRRFIQSN NLKAHILTHA NTNKNEQEGK
Source	Mammalian cell
Target Names	ZFP42
Protein Names	Recommended name: Zinc finger protein 42 homolog Short name= Zfp-42 Alternative name(s): Reduced expression protein 1 Short name= REX-1 Short name= hREX-1 Zinc finger protein 754
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.