



Recombinant Human Uncharacterized protein C20orf112 (C20orf112)

Product Code	CSB-MP856989HU
Abbreviation	C20orf112
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MY1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSDSTWMSAD PHLASSLSPS QDERMRSPQN LHSQEDDDSS SESGSGNGSS TLNPSTSSST QGDPAFPEMN GNGAVAPMDF TTAEDQPIN LCDKLPPATA LGTASYPSDG CGADGLRSRV KYGVKTTPEP PPYSSGSYDS IKTEVSGCPE DLTVGRAPTA DDDDDDHDDH EDNDKMNDSE GMDPERLKAF NMFVRLFVDE NLDRMVPISK QPKEKIQAI ESCSRQFPEF QERARKRIRT YLKSCRMRMK NGMEMTRPTP PHLTSAMAEN ILAAACESET RKAARKMRLE IYQSSQDEPI ALDKQHSRDS AAITHSTYSL PASSYSQDPV YANGGLNYSY RGYGALSSNL QPPASLQTGN HSNGPTDLSM KGGASTTSTT PTPTPSSTST SRPVPTAQLS PTEISAVRQL IAGYRESAAF LLRSADELEN LILQQN
Source	Mammalian cell
Target Names	NOL4L
Protein Names	Recommended name: Uncharacterized protein C20orf112
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.