



# Recombinant Human RUN domain-containing protein 3B (RUNDC3B)

<b>Product Code</b>	CSB-EP856993HU
<b>Abbreviation</b>	RUNDC3B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96NLO
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASRSLGGLS GIRGGGGGGG KKSLSARNAA VERRNLITVC RFSVKTLIDR SCFETIDDSS PEFNFAAIL EQILSHRLKE ISQSCRWLAH LQIPLQGQVT WFGYESPRSF WDYIRVACRK VSQNCICSIE NMENVSSRA KGRAWIRVAL MEKHLSEYIS TALRDFKTR RFYEDGAIVL GEEANMLAGM LLGLNAIDFS FCLKGEGLDG SFPVIDYTP YLKYIQSSDS ISSDEEELRT LGSSGSESS PENVGPPFLM DENSWFNKCK RVKQKYQLTL EQKGYLEELL RLRENQLSES VSQNKILLQR IEDSDLAHKL EKEQLEYIIV ELQDQLTVLK NNDLRSRQEL TAHLTNQWPS PGALDVNAVA LDTLLYRKHN KQWYEKSYQS LDQLSAEVSL SQTSLDPGQS QEGDGKQDTL NVMSEGKEDT PSELLGLCGSL TSVASYKSLT SLKSNDYLAS PTTEMTSPGL TPS
<b>Source</b>	E.coli
<b>Target Names</b>	RUNDC3B
<b>Protein Names</b>	Recommended name: RUN domain-containing protein 3B Alternative name(s): Rap2-binding protein 9 Rap2-interacting protein 9 Short name= RPIP-9
<b>Expression Region</b>	1-473
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.