



Recombinant Human RUN domain-containing protein 3B (RUNDC3B)

Product Code	CSB-MP856993HU
Abbreviation	RUNDC3B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NLO
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASRSLGGLS GIRGGGGGGG KKSLSARNAA VERRNLITVC RFSVKTLIDR SCFETIDDSS PEFNNFAAIL EQILSHRLKE ISQSCRWLAH LQIPLQGQVT WFGYESPRSF WDYIRVACRK VSQNCICSIE NMENVSSSRA KGRAWIRVAL MEKHLSEYIS TALRDFKTTR RFYEDGAIVL GEEANMLAGM LLGLNAIDFS FCLKGEGLDG SFPVIDYTP YLKYIQSSDS ISSDEEELRT LGSSGSESST PENVGPPFLM DENSWFNKCK RVKQKYQLTL EQKGYLEELL RLRENQLSES VSQNKILLQR IEDSDLAHKL EKEQLEYIIV ELQDQLTVLK NNDLRSRQEL TAHLTNQWPS PGALDVNAVA LDTLLYRKHN KQWYEKSYQS LDQLSAEVSL SQTSLDPGQS QEGDGKQDTL NVMSEGKEDT PSLGLCGSL TSVASYKSLT SLKSNDYLAS PTTEMTSPGL TPS
Source	Mammalian cell
Target Names	RUNDC3B
Protein Names	Recommended name: RUN domain-containing protein 3B Alternative name(s): Rap2-binding protein 9 Rap2-interacting protein 9 Short name= RPIP-9
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.