



Recombinant Human Serpin B11 (SERPINB11)

Product Code	CSB-EP856997HU
Abbreviation	SERPINB11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96P15
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGSLSTANVE FCLDVFKELN SNNIGDNIFF SLSLLYALS MVLLGARGET AEQLEKVLHF SHTVDSLKPG FKDSPKCSQA GRIHSEFGVE FSQINQPDSN CTLSIANRLY GTKTMAFHQQ YLSCSEKWYQ ARLQTVDFEQ STEETRKMIN AWVENKTNGK VANLFGKSTI DPSSVMVLVN IYFKGQRQN KFQVRETVKS PFQLSEGKNV TVEMMYQIGT FKLA FVKEPQ MQVLELPYVN NKLSMIILLP VGIANLKQIE KQLNSGTFHE WTSSSNMMER EVEVHLPRFK LEIKYELNSL LKPLGVTDLF NQVKADLSGM SPTKGLYLSK AIHKS YLDVS EEGTEAAAAT GDSIAVKSLP MRAQFKANHP FLFFIRHTHT NTILFCGKLA SP
Source	E.coli
Target Names	SERPINB11
Protein Names	Recommended name: Serpin B11
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.