



# Recombinant Human BRCA1-A complex subunit RAP80 (UIMC1)

<b>Product Code</b>	CSB-BP857019HU
<b>Abbreviation</b>	UIMC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96RL1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPRRKKKVKE VSESRNLEKK DVETTSSVSV KRKRRLLEDAF IVISDSDGEE PKEENGLQKT KTKQSNRAKC LAKRKIAQMT EEEQFALALK MSEQEAREVN SQEEEEEEELL RKAIAESLNS CRPSDASATR SRPLATGPSS QSHQEKT TDS GLTEGIWQLV PPSLFKGS HI SQGNEAEERE EPWDHTEKTE EEPVSGSSGS WDQSSQP VFE NVNVKSFDR C TGHS AEHTQC GKPQESTGRG SAFLKAVQGS GDTSRHCLPT LADAKGLQDT GGTVNYFWGI PFCPDGVDPN QYTKVILCQL EVYQKSLKMA QRQLLNKKG F GEPVLP RPPS LIQNECGQGE QASEKNECIS EDMGDEDKEE RQESRASDWH SKTKDFQESS IKSLKEKLLL EEPTTSHGQ SSQGIVEETS EEGNSVPASQ SVAALTSKRS LVLMPESAE EITVCPETQL SSSETFDLER EVSPGSRDIL DGVRIIMADK EVGNKEDAEK EVAISTFSS NQVSCPLCDQ CFPPTKIERH AMYCNGLMEE DTVLTRRQKE AKTKSDSGTA AQTSLDIDKN EKCYLCKSLV PFREYQCHVD SCLQLAKADQ GDGPEGSGRA CSTVEGKWQQ RLKNPKEKGH SEGRLLSFLE QSEHKTS DAD IKSSETGAFR VPSPGMEEAG CSREMQSSFT RRD LNESPVK SFVSI SEATD CLVDFKKQVT VQPGSRTRTK AGRGRRRK F
<b>Source</b>	Baculovirus
<b>Target Names</b>	UIMC1
<b>Protein Names</b>	Recommended name: BRCA1-A complex subunit RAP80 Alternative name(s): Receptor-associated protein 80 Retinoid X receptor-interacting protein 110 Ubiquitin interaction motif-containing protein 1
<b>Expression Region</b>	1-719
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.