



# Recombinant Clostridium acetobutylicum Glutamyl-tRNA (Gln) amidotransferase subunit A 2 (gatA2)

<b>Product Code</b>	CSB-BP857143DUB
<b>Abbreviation</b>	gatA2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q97EX8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium acetobutylicum (strain ATCC 824 / DSM 792 / JCM 1419 / LMG 5710 / VKM B-1787)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEILDMTVEQ LRNAILDKHL KSEDIVKAYF DNIKRNEPEI NAYITLCEDY ALKEAKDVVK KIANGDKVGR LAGIPIAIKD NICTDGIKTT CASKMLYDFV PPYDATVIKK LKAEDAIIG KVNMDDEFAMG SSTENSAFKI TKNPRDITRV PGGSSGGSAA VVAAKMAPIS LGSDTGGSIR QPAAFCGVVG LKPTYGLVSR FGLIAFASSL DQIGPLGKT V KDCAELLEVI SGEDELNTS SKKHEKEDYL EGIDDIGKGM KIGMPKEFLN DGLDPEIRKC IDDTIEKLKS LGAEVCEMSL PITEEGLSAY YIISSAEASS NLARFDGIRY GYRPDDFEDV YDLMETSRSR AFGDEVKRRRI MLGTYALSSG YYDAYYKRAL KLKKKIKNEF KEAFENYDLI LSPVSPVLPF KIGEKKADPL QMYLADIYTV NINLAGIPAI SLPCSVSKEG LPIGLQLLGP HFGEKKIFRA ARALEKER
<b>Source</b>	Baculovirus
<b>Target Names</b>	gatA2
<b>Protein Names</b>	Recommended name: Glutamyl-tRNA(Gln) amidotransferase subunit A 2 Short name= Glu-ADT subunit A 2 EC= 6.3.5.-
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.