



Recombinant Clostridium acetobutylicum Phosphatidylserine decarboxylase proenzyme 2 (psd2)

Product Code	CSB-EP857175DUB
Abbreviation	psd2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q97KW7
Product Type	Recombinant Protein
Immunogen Species	Clostridium acetobutylicum (strain ATCC 824 / DSM 792 / JCM 1419 / LMG 5710 / VKM B-1787)
Purity	>85% (SDS-PAGE)
Sequence	MSITIFNRDL KEIYYEKQYK NGQLKFLYNT ILGRMLLKLF ISTRLFSKIN AIFNNKGS I KKIEPFINEY KIDMSEYEK EYTSFDDFFT RKILEGKR SF SKEKSHLISP ADSKLMVYEI DDDLKMN IKN SIYTVGELLN DEKLSREYKN GTCLIFRLTV DDYHRYCFID DGSLKYRKVI NGR LHTVGPI SSKRYKVYSE NNREYSVLKT RNFGKVIQIE VGALLVGKIK NHSIKVFKKG DEKGYFCFGG
Source	E.coli
Target Names	psd2
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme 2 EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2 2. Phosphatidylserine decarboxylase beta chain 2
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.