



Recombinant Danio rerio Homeobox protein Hox-A1a (hoxa1a)

Product Code	CSB-MP857418DIL
Abbreviation	hoxa1a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q98SI1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MSTFLDFSSI SGGGDGGSGG SCSVRAFHGD HGLSTFQSSC AVRLNNSCSGD ERFMSNISSQ DVINSQPQQA GSYQSPGTLT ITYSAHPSYG TQSFCTGYNH YALNQDVESS VSFPQCGPLV YSGNISSTVV QHRHHRHGYS SGNVHLHGQF QYGSATYGNS SDQANLTFVA GCSNPLSPLH VPHHDACCSP LSDGVPTGQT FDWMKVKRNP PKTGKAGEYG FGGQPNTVRT NFSTKQLTEL EKEFHFNKYL TRARRVEIAA SLQLNETQVK IWFQNRMRKQ KKREKEGLLP KSLSEQKDGL EKTEDASEKS PSAPSTSPSPS PTVEAYSSN
Source	Mammalian cell
Target Names	hoxa1a
Protein Names	Recommended name: Homeobox protein Hox-A1a Short name= Hox-A1
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.