



Recombinant Human ETS-related transcription factor Elf-4 (ELF4)

Product Code	CSB-YP857453HU
Abbreviation	ELF4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99607
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAITLQPSDL IFEFASNGMD DDIHQLEDPS VFPAVIVEQV PYPDLLHLYS GLELDDVHNG IITDGTLCMT QDQILEGSFL LTDDNEATSH TMSTAEVLLN MESPSDILDE KQIFSTSEML PDSDPAPA VT LPNYLFPASE PDALNRAGDT SDQEGHSLEE KASREESAKK TGKSKKRIRK TKGNRSTSPV TDPSIPIRKK SKDGGKSTIY LWEFLLALLQ DRNTCPKYIK WTQREKGIFK LVDSKAVSKL WGKQKNK PDM NYETMGRALR YYYQRGILAK VEGQRLVYQF KEMPKDLVVI EDEDESSEAT AAPPQASTAS VASASTTRRT SSRVSSRSAP QGKGSSSWEK PKIQHVGLQP SASLELGPSL DEEIPTTSTM LVSPAEGQVK LTKAVSASSV PSNIHLGVAP VGSGSALTLQ TIPLTTVL TN GPPASTTAPT QLVLQSV PAA STFKDTFTLQ ASFPLNASFQ DSQVAAPGAP LILSGLPQLL AGANRPTNPA PPTVTGAGPA GPSSQPPGTV IAAFIRTS GT TAAPRVKEGP LRSSSYVQGM VTGAPMEGLL VPEETLRELL RDQAHLQPLP TQVVSRGSHN PSELLGNQTLS PPSRPTVGLT PVAEELSSG SGSLLEMAEPS VTTSGSLLTR SPTPAPFSPF NPTSLIKMEP HDI
Source	Yeast
Target Names	ELF4
Protein Names	Recommended name: ETS-related transcription factor Elf-4 Alternative name(s): E74-like factor 4 Myeloid Elf-1-like factor
Expression Region	1-663
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.