



# Recombinant Mouse Forkhead box protein P3 (Foxp3)

<b>Product Code</b>	CSB-EP857473MO-B
<b>Abbreviation</b>	Foxp3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99JB6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPNPRPAKPMAPSLALGPSPGVLPWKTAPKGSELLGTRGSGGPFQGRDLR SGAHTSSSLNPLPPSQLQLPTVPLVMVAPSGARLGPSPHLQALLQDRPHFMH QLSTVDAHAQTPVLQVRPLDNPAMISLPPPSAATGVFSLKARPGLPPGINVASL EWVSREBALLCTFPRSGTPRKDSNLLAAPQGSYPLLANGVCKWPGCEKVFE PEEFLKHCQADHLLDEKGAQCLLQREVVSLEQQLELEKEKLGAMQAHLA KMALAKAPSVASMDKSSCCIVATSTQGSVLPAWSAPREAPDGGLFAVRRHLW GSHGNSSFPEFFHNMDYFKYHNMRPPFTYATLIRWAILEAPERQRTLNEIYHW FTRMFAYFRNHPATWKNAIRHNLHLKCFVRVESEKGAVWTVDEFEFRKKR
<b>Source</b>	E.coli
<b>Target Names</b>	Foxp3
<b>Protein Names</b>	Recommended name: Forkhead box protein P3 Alternative name(s): Scurfin
<b>Expression Region</b>	1-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.