



Recombinant Mouse Pre-mRNA-processing factor 19 (Prpf19)

Product Code	CSB-EP857484MO
Abbreviation	Prpf19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99KP6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SLICISISNE VPEHPCVSPV SNHVYERRLI EKZIAENGTD PINNQPLSEE QLIDIKVAHP IRPKPPSATS IPAILKALQD EWDAVMLHSF TLRQQLQTRR QELSHALYQH DAACRVIARL TKEVTAAREA LATLKPQAGL IVPQAVPSSQ PSVVGAGEPM DLGELVGMTP EIIQKLQDKA TVLTTERKKR GKTVPEELVK PEELSKYRQV ASHVGLHSAS IPGILALDLC PSDTNKILTG GADKNVVVFD KSTEQILATL KGHTKKVTSV VFHPSQELVF SASPDATIRI WSPNTSCVQ VVRAHESAVT GLSLHATGDY LLSSDDQYW AFSDIQTGRV LTKVTDETSG CSLTCAQFHP DGLIFGTGTM DSQIKIWDLK ERTNVANFPG HSGPITSIAF SENGYYLATA ADDSSVKLWD LRKLKNFKTL QLDNNFEVKS LIFDQSGTYL ALGGTDVQIY ICKQWTEILH FTEHSGLTTG VAFGHHAKFI ASTGMDRSLK FYSL
Source	E.coli
Target Names	Prpf19
Protein Names	Recommended name: Pre-mRNA-processing factor 19 Alternative name(s): Nuclear matrix protein 200 PRP19/PSO4 homolog Senescence evasion factor
Expression Region	2-504
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.