



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-L1-1 (CYCL1-1)

Product Code	CSB-MP857683OFG
Abbreviation	CYCL1-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9AS36
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MIYTAIDTFY LTDEQLRDSP SRKDGIDEAT ETALRVYGCD LIQESGILLK LPQAVMATAQ VLFHRFYCKK SFVRFSVKRV AASCVWLAGK LEESPRRSKH IIIVFHRMEC RRENVPIEHL DVFSKKYSDL KHDLVRTERH LLKEMGFICH VEHPHKFISN YLATLEAPEL TQEAWNLANL SLRTTLCVRF KSEVVACGVV YAAARRHGVP LPEDPPWWNV FDADEAGIQE VCRVLAHLYS LPKSQYIQVY KDNDSTFTHRR TSDTNASKES PATTVASDKG TPVPSSSSQE KDALIKAGSD KVKEKGDDDG KTLPSEPNGK EGPAVNLKSE KSESNVDRSR ERERDRSRGR DRDSRGRDSD RDSKGRDSDR ERERDREADR DRQRRHHSKD RSSGYSDKEK SRHRSSDRG DHYSSHSSRD KDRHRRQ
Source	Mammalian cell
Target Names	CYCL1-1
Protein Names	Recommended name: Cyclin-L1-1 Short name= Cycl1;1
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.