



Recombinant Mouse Ubiquitin carboxyl-terminal hydrolase BAP1 (Bap1)

Product Code	CSB-MP858174MO
Abbreviation	Bap1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99PU7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNKGWLELES DPGLFTLLVE DFGVKGQVE EIYDLQSKCQ GPVYGFIFLF KWIEERRSRK KVSTLVDDTS VIDDDIVNNM FFAHQLIPNS CATHALLSVL LNCNSVDLGP TLSRMKDFTK GFSPESKGYA IGNAPELAKA HNSHARPEPR HLPEKQNGLS AVRTMEAFHF VSYVPITGRL FELDGLKVYP IDHGPWGEDE EWTDKARRVI MERIGLATAG EPYHDIRFNL MAVVPDRRIK YETRLHVLKV NRQTVLEALQ QLIRVTQPEL IQTHKSQESQ LPEESKPASS KSPLGLEAGR TPVASECTQT DGAEVAVGSC PQTTHSPSS KCKLVVKKPPG SSLNGVPPNP APIVQRLPAF LDNHNYAKSP MQEEEDLAAG VGRSRVPVRA PQQYSEDEDD YEDEDVQN TNPAIRYKRK GTGKPGSLSN SSDGQLSVLQ PNTINVLTEK LQESQKDLV PLSIKTSSGA GSPAVAVPTH SQPSPTPSNE STDTASEIGS AFNSPLRSPI RSANPTRPSS PVTSHISKVL FGEDDSLLRV DCIRYNRAVR DLGPVISTGL LHLAEDGVLS PLALTEGGKG SSPSTRSSQG SQGSSGLEEK EVVEVTESRD KPGLNRSSEP LSGEKYSPE LLALLKVEA EIANYEACLK EEVEKRKKFK IDDQRRTHNY DEFICTFISM LAQEGMLANL VEQNISVRRR QGVSIGRLHK QRKPDRRKR RPYKAKRQ
Source	Mammalian cell
Target Names	Bap1
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase BAP1 EC= 3.4.19.12 Alternative name(s): BRCA1-associated protein 1 Ubiquitin C-terminal hydrolase X4 Short name= UCH-X4
Expression Region	1-728
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.