



# Recombinant *Oryza sativa* subsp. *japonica* Ribosome production factor 2 homolog (Os01g0513800, LOC\_Os01g33030)

<b>Product Code</b>	CSB-BP858343OFG
<b>Abbreviation</b>	Os01g0513800, LOC_Os01g33030
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9AWM9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVAAIRVPRS QRAKRELLKH APKLVETGKK TLILHGTKTS AVLNSVLADL FHLKRDNAVK YSKKNDNIRP FESGGETSLE FFSLKTDCSL IVYGSHSKKR PNNLVLGRTY DHHIYDLVEV GVENYKSIES YVYDKKLAPK LGSKPFFAFI GEHFESVEEL KHLKEVLLDL FKGEVVENLN LAGVDRVFVC TAISPTTVYM MHCALRLKRS GTSIPRIELV EVGPSMDLVV RRHRYPVESL KKEAMKTADH AKKMKNVTKD PVHGKLGKVY IPDQQIAKMS LSNDVKGLKR ERREAKKNKD HSKKQKINPE
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os01g0513800
<b>Protein Names</b>	Recommended name: Ribosome production factor 2 homolog Alternative name(s): Brix domain-containing protein 1 homolog Ribosome biogenesis protein RPF2 homolog
<b>Expression Region</b>	1-310
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.