



# Recombinant Human Mannose-6-phosphate isomerase (MPI)

<b>Product Code</b>	CSB-EP014754HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P34949
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AAPRVFPLS CAVQQYAWGK MGSNSEVARL LASSDPLAQI AEDKPYAELW MGTHPRGDAK ILDNRSQKT LSQWIAENQD SLGSKVKDTF NGNLPFLFKV LSVETPLSIQ AHPNKELAEK LHLQAPQHYP DANHKPEMAI ALTPFQGLCG FRPVEEIVTF LKKVPEFQFL IGDEAATHLK QTMSHDSQAV ASSLQSCFSH LMKSEKKVVV EQLNLLVKRI SQQAAAGNNM EDIFGELLQ LHQQYPGDIG CFAIYFLNLL TLKPGEAMFL EANVPHAYLK GDCVECMACS DNTVRAGLTP KFIDVPTLCE MLSYTPSSSK DRLFLPTRSQ EDPYLSIYDP PVPDFTIMKT EVPGSVTEYK VLALDSASIL LMVQGTVIAS TPTTQTPIPL QRGGVLFIGA NESVSLKLTE PKDLLIFRAC CLL
<b>Source</b>	E.coli
<b>Target Names</b>	MPI
<b>Protein Names</b>	Recommended name: Mannose-6-phosphate isomerase EC= 5.3.1.8 Alternative name(s): Phosphohexomutase Phosphomannose isomerase Short name= PMI
<b>Expression Region</b>	2-423
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Phosphomannose isomerase catalyzes the interconversion of fructose-6-phosphate and mannose-6-phosphate and plays a critical role in maintaining the supply of D-mannose derivatives, which are required for most glycosylation reactions. Mutations in the MPI gene were found in patients with carbohydrate-deficient glycoprotein syndrome, type Ib.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.