



Recombinant *Saccharomyces cerevisiae* Pre-rRNA-processing protein PNO1 (PNO1)

Product Code	CSB-EP858700SVG-B
Abbreviation	PNO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99216
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVAPTALKKA TVTPVSGQDG GSSRIIGINN TESIDEDDDD DVLLDDSDNN TAKEEVEGEE GSRKTHESKT VVDDQ GKPR FTSASKTQGN KIKFESRKIM VPPHRMTPLR NSWTKIYPPL VEHLKLQVRM NLKTKSVELR TNPKFTTDPG ALQKGADFIK AFTLGFDLDD SIALLRLLDD YIETFEVKDV KTLTGDHLRSR AIGRIAGKDG KTKFAIENAT RTRIVLADSK IHILGGFTHI RMARES VVSL ILGSPPGKVY GNLRTVASRL KERY
Source	E.coli
Target Names	PNO1
Protein Names	Recommended name: Pre-rRNA-processing protein PNO1 Alternative name(s): Partner of NOB1 Ribosomal RNA-processing protein 20
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.