



Recombinant *Saccharomyces cerevisiae* Trafficking protein particle complex subunit 33 (TRS33)

Product Code	CSB-BP858709SVG
Abbreviation	TRS33
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99394
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SSTHSNNVG HPQSSPQGGL TEQQRAQQQY QIFENSLPKV SQSVYQMLLN EMVPLAMGIE RQISGDVISS DSNVTSENGN INNMIKRLKI EEHHTVDIIR SHNLIHELYK ADEEKEKVL ARLRNIGFQI GLKLSSELLIF SNNPNLKFKE MDLLLIMKFI CRDVWKQIFG KQIDNLKTNH RGTFYLLDYD YRPIQSFSLE EDAKNEELKM IEPFLEIPVG IIRGVLSSLG YSSEEVICLA SFIDRPTDRP KTAFPKGVSF HVQVTMPQ
Source	Baculovirus
Target Names	TRS33
Protein Names	Recommended name: Trafficking protein particle complex subunit 33 Short name= TRAPP subunit 33 Alternative name(s): Transport protein particle 33 kDa subunit
Expression Region	2-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.