



Recombinant Mouse ADP-ribosylation factor GTPase-activating protein 2 (Arfgap2)

Product Code	CSB-YP858747MO
Abbreviation	Arfgap2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99K28
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AASPSKTEI QTIFKRLRAI PTNKACFDCG AKSPSWASIT YGVFLCIDCS GVHRSLGVHL SFIRSTELDS NWSWLQLRCM QVGGNANATA FFRQHGCMA DANTKYTSRA AQMYREKIRQ LGSAALTRHG TDLWIDSMNS APSHSPEKKD SDFTEHTQA PAWDTAATDP SGTQQPALPS ESSSLAQPEQ GPNTDLLGTS PQASLELKSS IIGKKKAAA KKGLGAKKGL GAQKVSNSQSF TEIERQAQVA EKLREQQAAD AKKQAEESMV ASMRLAYQEL QIDRKKEEKK LQNLEGKKRE QAERLGMGLV SRSSISHSVL SEMQMIEQET PLSAKSSRSQ LDFDDVGTG ASGPPKYKDN PFSLGETFGS RWDSDAAWGM DRVEEKEPEV TISSIRPISE RTASRREVET RSSGLESSEA RQKFAGAKAI SSDMFFGREV DSEYEARSRL QQLSGSSAIS SSDLFGNMDG AHGGGTVSLG NVLPTADIAQ FKQGVKSVAG KMAVLANGVM NSLQDRYGSY
Source	Yeast
Target Names	Arfgap2
Protein Names	Recommended name: ADP-ribosylation factor GTPase-activating protein 2 Short name= ARF GAP 2 Alternative name(s): GTPase-activating protein ZNF289 Zinc finger protein 289
Expression Region	2-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.