



Recombinant Transcriptional regulatory protein ZraR (zraR)

Product Code	CSB-MP858944KBE
Abbreviation	zraR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9APD9
Product Type	Recombinant Protein
Immunogen Species	Klebsiella oxytoca
Purity	>85% (SDS-PAGE)
Sequence	MSGEQVDILV VDDDISHCTI LQALLRGWGY RVALANGLQ ALEKVREKVF DLVLCDIRMA EMDGIETLKE IKTFNPSIPV LIMTAYSSVD TAVEALKSGA LDYLIKPLDF DKLQLTLSEA LAHTRLSESP VTETPAAQFG MVGDSPAMRA LLNNITLVAP SDATVLIHGE SGTGKELVAR ALHASSARSR RPLVILNCAA LNESLLESEL FGHEKGAFTG ADKRREGRFV EADGGTLFLD EIGDISPLMQ VRLLRAIQER EVQRVGSNQT LSV DVRLIAA THRD LAEEVS AGRFRQDLYY RLNVVTIDMP PLRHRREDIP PLARYFLQRY AERNRKAVQG FTPQAMDLLI HYAWPGNIRE LENAVERAVV LLTGEYISER ELPLAITGTP VADAPHGDDS IQPLVEVEKE AILAAALERTG GNKTEAARRL GITRKTLLAK LSR
Source	Mammalian cell
Target Names	zraR
Protein Names	Recommended name: Transcriptional regulatory protein ZraR
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.