



# Recombinant *Oryza sativa* subsp. japonica AP2/ERF and B3 domain-containing protein Os01g0141000 (Os01g0141000, LOC\_Os01g04800)

<b>Product Code</b>	CSB-EP858960OFG-B
<b>Abbreviation</b>	Os01g0141000, LOC_Os01g04800
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9AWS0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGVVSFSSTS SGASTATTES GGAVRMSPEP VVAVAAAAQQ LPVVKGVDSA DEVVTSRPAA AAAQSSRYK GVPQPNGRW GAQIYERHAR VWLGTFPDEE AARAYDVAA LRYRGRDAAT NFPGAAASAA ELAFLAAHSK AEIVDMLRKH TYADELRQGL RRGRGMGARA QPTPSWAREP LFEKAVTPSD VGKLNRLVVP KQHAEKHFPL RRASSDSAS AAATGKGVLL NFEDGEKVV RFRYSYWNSS QSYVLTKGWS RFVREKGLRA GDTIVFSRSA YGPDKLLFID CKKNNAAT TTCAGDERPT TSGAEPVVR LFGVDIAGGD CRKRERAVEM GQEVFLLKRQ CVVHQRTPAL GALLL
<b>Source</b>	E.coli
<b>Target Names</b>	Os01g0141000
<b>Protein Names</b>	Recommended name: AP2/ERF and B3 domain-containing protein Os01g0141000
<b>Expression Region</b>	1-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.