



Recombinant *Oryza sativa* subsp. *japonica* AP2/ERF and B3 domain-containing protein Os01g0141000 (Os01g0141000, LOC_Os01g04800)

Product Code	CSB-MP858960OFG
Abbreviation	Os01g0141000, LOC_Os01g04800
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9AWS0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGVVSFSSTS SGASTATTES GGAVRMSPEP VVAVAAAAQQ LPVVKGVDSA DEVVTSRPAA AAAQSSRYK GVPQPNGRW GAQIYERHAR VWLGTFPDEE AARAYDVAA LRYRGRDAAT NFPGAAASAA ELAFLAAHSK AEIVDMLRKH TYADELRQGL RRGRGMGARA QPTPSWAREP LFEKAVTPSD VGKLNRLVVP KQHAEKHFPL RRASSDSAS AAATGKGVLL NFEDGEKVVW RFRYSYWNSS QSYVLTKGWS RFVREKGLRA GDTIVFSRSA YGPDKLLFID CKKNNAAAAT TTCAGDERPT TSGAEPRVVR LFGVDIAGGD CRKRERAVEM GQEVFLLKRQ CVVHQRTPAL GALLL
Source	Mammalian cell
Target Names	Os01g0141000
Protein Names	Recommended name: AP2/ERF and B3 domain-containing protein Os01g0141000
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.