



# Recombinant *Saccharomyces cerevisiae* Transposon Ty3-I Gag polyprotein (TY3A-I)

<b>Product Code</b>	CSB-MP859088SVG
<b>Abbreviation</b>	TY3A-I
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99219
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SFMDQIPGG GNYPKLPVEC LPNFPIQP SL TFRGRNDSHK LKNFISEIML NMSMISWPND ASRIVYCRRH LLNPAAQWAN DFVQEQGILE ITFDTFIQGL YQHFYKPPDI NKIFNAITQL SEAKLGIERL NQRFRKIWDR MPPDFMTEKA AIMTYTRLLT KETYNIVRMH KPETLKDAME EAYQTTALTE RFFPGFELDA DGDTIIGATT HLQEEYDSY DSEDNLTQNR YVHTVTRRRS YNKPMSNHRN RRNNNASREE CIKNRLCFYC KKEGHLNEC RARKASSNRS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TY3A-I
<b>Protein Names</b>	Recommended name: Transposon Ty3-I Gag polyprotein Alternative name(s): Gag3 Transposon Ty3-2 protein A Short name= TY3A Cleaved into the following 3 chains: 1. Capsid protein Short name= 2. CA Alternative name(s): p24
<b>Expression Region</b>	2-290
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.