



Recombinant *Saccharomyces cerevisiae* Family of serine hydrolases 3 (FSH3)

Product Code	CSB-MP859095SVG
Abbreviation	FSH3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99369
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSEKKKVLML HGFVQSDKIF SAKTGGLRKN LKKLGYDLYY PCAPHSIDKK ALFQSESEKG RDAAKEFNST ATSDEVYGWF FRNPESFNSF QIDQKVFNYL RNYVLENGPF DGVIGFSQGA GLGGYLVTFD NRILNLTDEQ QPALKFFISF SGFKLEDQSY QKEYHRIVQ PSLHVRGELD EVVAESRIMA LYESWPDNKR TLLVHPGAHF VPNSKPFVSQ VCNWIIQGITS KEGQEHNAPQ EVDRKQFDKP QLEDDLLDMI DSLGKL
Source	Mammalian cell
Target Names	FSH3
Protein Names	Recommended name: Family of serine hydrolases 3 EC= 3.1.-.-
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.