



Recombinant Human Protein kinase C-binding protein NELL2 (NELL2)

Product Code	CSB-BP859098HU
Abbreviation	NELL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99435
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> LGVDPSLQI DVLTELELGE STTGVRQVPG LHNGTKAFLF QDTPRSIKAS TATAEQFFQK LRNKHEFTIL VTLKQTHLNS GVILSIHHL D HRYLELESSG HRNEVRLHYR SGSHRPHTEV FPYILADDKW HKLSLAISAS HLILHIDCNK IYERVVEKPS TDLPLGTTFW LGQRNNAHGY FKGIMQDVQL LVMPQGFI AQ CPDLNRTCPT CNDFHGLVQK IMELQDILAK TSAKLSRAEQ RMNRLDQCYC ERTCTMKGTT YREFESWIDG CKNCTCLNGT IQCETLICPN PDCPLKSALA YVDGKCKKEC KSICQFQGR T YFEGERNVY SSSGVCVLYE CKDQTMKLVE SSGCPALDCP ESHQITLSHS CCKVCKGYDF CSERHNCMEN SICRNLDRA VCSCRDGFRA LREDNAYCED IDECAEGRHY CRENTMCVNT PGSFMCICKT GYIRIDDYSC TEHDECITNQ HNCDENALCF NTVGGHNCVC KPGYTGNGTT CKAFCKDGCR NGGACIAANV CACPQGFTGP SCETDIDECS DGFVQCDSRA NCINLPGWYH CECRDGYHDN GMFSPSGESC EDIDECGTGR HSCANDTICF NLDGGYDCRC PHGKNCTGDC IHDGKVKHNG QIWLLENDRC SVCSCQNGFV MCRRMVCDC E NPTVDLFCCP ECDPRLSSQC LHQNGETLYN SGDTWVQNCQ QCRCLQGEVD CWPLPCPDVE CEFSILPENE CCPRCVTDPC QADTIRNDIT KTCLDEMNVV RFTGSSWIKH GTECTLCQCK NGHICCSVDP QCLQEL </p>
Source	Baculovirus
Target Names	NELL2
Protein Names	Recommended name: Protein kinase C-binding protein NELL2 Alternative name(s): NEL-like protein 2 Nel-related protein 2
Expression Region	22-816
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a glycoprotein containing several von Willebrand factor C



domains and epidermal growth factor (EGF)-like domains. The encoded protein acts as a homotrimer and is found in the cytoplasm. Several variants encoding a few different isoforms exist, and at least one isoform appears to be a secreted protein. Studies in mouse suggest that this protein plays a role in neural cell growth and differentiation as well as in oncogenesis.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.