



Recombinant Human Cell growth regulator with RING finger domain protein 1 (CGRRF1)

Product Code	CSB-BP859107HU
Abbreviation	CGRRF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99675
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAVFLVTLY EYSPLFYIAV VFTCFIVTTG LVLGWFGWDV PVILRNSEET QFSTRVFKKQ MRQVKNPFG L EITNPSSASI TTGITLTTDC LEDSLLT CYW GCSVQKLYEA LQKHVYCFRI STPQALEDAL YSELYQEY FIKKDSKEEI YCQLPRDTKI EDFGTVPRSR YPLVALLTLA DEDDREIYDI ISMVSVIHIP DRTYKLSCRI LYQYLLLAQG QFHDLKQLFM SANNNFTPSN NSSSEEKNTD RSLLEKVGLS ESEVEPSEEN SKDCVVCQNG TVNWWLLPCR HTCLCDGCVK YFQQCPMCRQ FVQESFALCS QKEQDKDKPK TL
Source	Baculovirus
Target Names	CGRRF1
Protein Names	Recommended name: Cell growth regulator with RING finger domain protein 1 Alternative name(s): Cell growth regulatory gene 19 protein RING finger protein 197
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.