



Recombinant Human Evolutionarily conserved signaling intermediate in Toll pathway, mitochondrial (ECSIT)

Product Code	CSB-EP859414HU-B
Abbreviation	ECSIT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQ95
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SS EQSLVSPPE PRQRPTKALV PFEDLFGQAP GGERDKASFL QTVQKFAEHS VRKRGHIDFI YLALRKMREY GVERDLAVYN QLLNIFPKEV FRPRNIIQRI FVHYPRQQEC GIAVLEQMEN HGVMFNKETE FLLIQIFGRK SYPMLKLVRL KLWFRPFMNV NPFPVPRDLP QDPVELAMFG LRHMEPDLA RVTIYQVPLP KDSTGAADPP QPHIVGIQSP DQQAALARHN PARPVFVEGP FSLWLRNKCV YYHILRADLL PPEEREVEET PEEWNLYYPM QLDLEYVRS WDNYEFDINE VEEGPVFAMC MAGAHDQATM AKWIIQGLQET NPTLAQIPVV FRLAGSTREL QTSSAGLEEP PLPEDHQEED DNLQRQQQGG S
Source	E.coli
Target Names	ECSIT
Protein Names	Recommended name: Evolutionarily conserved signaling intermediate in Toll pathway, mitochondrial Alternative name(s): Protein SITPEC
Expression Region	49-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.