



Recombinant Pelodiscus sinensis japonicus L-lactate dehydrogenase B chain (LDHB)

Product Code	CSB-EP859481PLAX
Abbreviation	LDHB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q98SK9
Product Type	Recombinant Protein
Immunogen Species	Pelodiscus sinensis japonicus (Chinese soft-shelled turtle)
Purity	>85% (SDS-PAGE)
Sequence	ATLQEKLIT PIVAGSTTPN NKITVVGVGQ VGMACAISIL GKGLCDELAL VDVLEDKLGK EMMDLQHGS LFLQTHKIVAD KDLRVTANSK IVVVTAGVRQ QEGESRLNLV QRNVNVFKFI IPQIMKYSNP CTILVVSNPV DILTYVTWKL SGLPKHRVIG SGCNLD SARF RHLMAEKLGI HPTSCHGWIL GEHGDSSVAV WSGVNVAGVS LQELNPAMGT DKDSENWKEV HKMVVDSAYE VIKLGKGTNW AIGLSVADLI ESYVKNLCRV HPVSTMVKGM YGIENEVFLS LPCVLSASGL TSVINQKLKD EEVAQLRKSA DTLWSIQKDL KDL
Source	E.coli
Target Names	LDHB
Protein Names	Recommended name: L-lactate dehydrogenase B chain Short name= LDH-B EC= 1.1.1.27
Expression Region	2-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.