



Recombinant Danio rerio Eukaryotic translation initiation factor 4E-binding protein 3-like (eif4ebp3l)

Product Code	CSB-EP859485DIL
Abbreviation	eif4ebp3l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q98TT6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSTNTQQSKS CPIPTRVLHL KDWSQLPDCY SQTPGGTLFS TTPGGTRIIY DRKFLLDCRN SPIARTPPCC LPQIPGVTIP SLHPVSKLQE LKEELEEEKE LAADDSQFEM DI
Source	E.coli
Target Names	eif4ebp3l
Protein Names	Recommended name: Eukaryotic translation initiation factor 4E-binding protein 3-like Short name= 4E-BP3 Short name= eIF4E-binding protein 3 Alternative name(s): Eukaryotic translation initiation factor 4E-1A-binding protein Short
Expression Region	1-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.