



Recombinant *Saccharomyces cerevisiae* Protein HSH49 (HSH49)

Product Code	CSB-EP859497SVG
Abbreviation	HSH49
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99181
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNYSADSGNT VYVGNIDPRI TKEQLYELFI QINPVLRIKY PKDKVLQAYQ GYAFIEFYNQ GDAQYAIIKIM NNTVRLYDRL IKVRQVTNST GTTNLPSNIS KDMILPIAKL FIKNLADSID SDQLVKIFNK FGKLIREPEI FYLSNGKLKC AYVYFEDFEK ADLAIKSLNN QLVANNRITV DYAFKENGKG NAKYGDDVDR LLNKEALKHN MLK
Source	E.coli
Target Names	HSH49
Protein Names	Recommended name: Protein HSH49
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.