



Recombinant Mouse Cytoplasmic protein NCK1 (Nck1)

Product Code	CSB-MP859559MO
Abbreviation	Nck1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99M51
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AEEVVVAK FDYVAQQEQE LDIKKNERLW LLDDSKSWWR VRNSMNKTGF VPSNYVERKN SARKASIVKN LKDTLGIGKV KRKPSVPTA SPADDSFVDP GERLYDLNMP AFVKFNMAE REDELSIKG TKVIVMEKCS DGWWRGSYNG QIGWFPSNYV TEEGDSPLGD HVGSLSEKLA AVVNNLNTGQ VLHVQALYP FSSNDEELN FEKGDVMDVI EKPENDPEWW KCRKINGMVG LVPKNYVTIM QNNPLTSGLE PAPPQCDYIR PSLTGKFAGN PWYYGKVTRH QAEMALNERG HEGDFLIRDS ESSPNDFSVS LKAQGKKNHF KVQLKETVYC IGQRKFSTME ELVEHYKKAP IFTSEQGEKL YLVKHLS
Source	Mammalian cell
Target Names	Nck1
Protein Names	Recommended name: Cytoplasmic protein NCK1 Alternative name(s): NCK adaptor protein 1 Short name= Nck-1
Expression Region	2-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.