



Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor BRR1 (BRR1)

Product Code	CSB-BP859917SVG
Abbreviation	BRR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99177
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKRGESQAPD AIFGQSRAFA LSDSSVNPDV IEYLKSVRQE ALRTNAISIK NHMNLQKRTR HKSSMYDDED EGALKRHAIS PSLIRLQRNV EIWVRWFNSV KATVLTNAYE FTGYEDETLD LLLLFLKNYL EDMPSKCTTV EKISVLNQH SFPEKAEKE ENLQIDEEWA KNILVRLEKT KIDSVEDVKK VITEGDKHEL VGYNQWFQYL INNEPQHHTF HEKITSKQLW VLIKYSNTW IKEIHKKGRH YRRLQDWLFY ILVHTPERVT AEYTSILRDL GKKCLELIQK KPVEAHENKI TLPKEMAELN VEIPAAVENM TITELTVSVI AVNYGQKDLI E
Source	Baculovirus
Target Names	BRR1
Protein Names	Recommended name: Pre-mRNA-splicing factor BRR1 Alternative name(s): Bad response to refrigeration protein 1
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.